



# caribou<sup>TM</sup>

The Caribou is our new standing frame. The frame supports a prone or supine standing position. With its ability to switch from a horizontal to a vertical position, the Caribou is designed to accommodate a wide range of users.



A good standing position can be achieved for people with, for instance, spasticity, contractions and circulatory deficiencies.

The Caribou is built with a base board and a series of optional upper body supports to provide the optimum freedom of movement. Further accessories are available to meet the postural support needs of the user.

The Caribou is available in four sizes and suitable for children and young people.



0° - 87°



max.100 kg



51-75 cm



75-111 cm



### Easy transfer

The Caribou tilts from vertical to horizontal. The horizontal position makes it easier for the carer and less stressful for the user to be placed into position. The Caribou is designed with a raised frame for easy use of floor hoist




### Anatomic design and soft mouldable PU

New and soft PU material for optimum comfort. The shape can be moulded to the body using additional wedges to shape the PU. The Caribou is designed with smooth edges for better safety and it is easy to clean



### Accessories are easily interchangeable for carers convenience



### Easy tilt

The Caribou is supplied with an easy tilt mechanism, which allows the carer to have both hands free




## Standard

The Caribou can be ordered as either manual or power version



## Frame colour

Available in different colours




Yellow for size 1



Red for size 2



Blue for size 3



Metallic grey for size 4


## Foot plates

Parted foot supports makes it possible to adjust individual foot heights and angles.




## Tilt mechanism

Goes from minus 3 degrees to horizontal position





Foot operated gas spring



with power actuator

Sizes 3 and 4 are recommended in the power versions to ease tilting



## Measurements

	Size 1	Size 2	Size 3	Size 4
Frame width (G)	51 cm (20")	57 cm (22½")	65 cm (25¼")	75 cm (29¾")
Frame length, total (H)	75 cm (29¾")	92 cm (36")	111 cm (43¼")	111 cm (43¾")
Height, vertical (I)	46-57 cm (18-22½")	53-64 cm (20¾-25¼")	68-79 cm (26¾-31")	88-99 cm (34¾-39")
Height, horizontal (I1)	70 cm (27½")	70 cm (27½")	80 cm (31½")	80 cm (31½")
Height from floor to frame (I2)	13 cm (5")	13 cm (5")	13 cm (5")	13 cm (5")
Basis plate width, top-bottom (J+J1)	25-32 cm (9¾-12½")	30-38 cm (11¾-15")	35-46 cm (13¾-18")	40-55 cm (15¾-21¾")
Basis plate height (K)	43 cm (17")	50 cm (19¾")	65 cm (25½")	85 cm (33½")
Angling	0 - 87° (0 - 87°)	0 - 87° (0 - 87°)	0 - 87° (0 - 87°)	0 - 87° (0 - 87°)
Foot plate angling	-10 to 10°	-10 to 10°	-10 to 10°	-10 to 10°
Weight w/o power	20.5 kg (45.2lb)	24 kg (53lb)	30.5 kg (67.2lb)	35 kg (77.2lb)
Weight with power	22.5 kg (49.6lb)	25.5kg (56.2lb)	32.5 kg (71.7lb)	37 kg (81.6 lb)
Max. user weight / load	40 kg (88.2lb)	60 kg (132.2lb)	80 kg (176.4lb)	100 kg (220.5lb)



## ACCESSORIES



### Top plate

861101-x

Size	Measurements
1	W1:15, W2:18, H:20 cm
2	W1:19, W2:23, H:25 cm
3	W1:22, W2:27, H:32 cm
4	W1:25, W2:31, H:40 cm



### Chest strap

86101-x

Size	Measurements
1+2	W: 42-63 cm
3	W: 48-77 cm
4	W: 58-96 cm



### Cross straps

86102-x The cross straps should be positioned so that the top belt is level with the sacrum

Size	Measurements
1	38-49 cm
2	43-55 cm
3	53-70 cm
4	63-86 cm



### Side supports, set

861102-x The cushion can be turned for either high or deep side support

Size	Measurements
1	H:14 x W:9½ cm
2	H:18 x W:13 cm
3	H:18 x W:13 cm



### Straps for side supports

86103-x

Hip support close to the body

Size	Measurements
1	26-32 cm
2	31-41 cm
3	34-50 cm



### Shoulder support, set

861126-x

Especially for supine standing

Size	Measurements
1	W:7½, H:19, D:5 cm
2	W:8½, H:24, D:5 cm
3	W:10, H:31½, D:5 cm
4	W:11, H:39, D:5 cm



### Wedges for side slots

86107-x

Mould the plate close to the body



### Head support plate

861116-x

For the user that needs light support of the head

Size	Measurements
1	W1:15, W2:18, H:20 cm
2	W1:19, W2:23, H:25 cm



### Head support

All top plates and extra cushions are prepared for head rest mounting

See the Accessories section at the back of the catalogue, for further information on options and measurements



### Extra cushions

861106-x

Max. 2 extra cushions can be mounted

Size	Measurements
1	W:25 x H:8 cm
2	W:30 x H:10 cm
3	W:35 x H:10 cm
4	W:40 x H:10 cm



### Back support

861115-x

Gives a firm support over sacrum when used as a prone stander

Two sizes



### Pommel

861103-x

Size 1 for supine standing  
Size 2 for prone standing

Size	Measurements
1	W1:5, W2:7, H:15, D:17 cm
2	W1:15, W2:7, H:15, D:9 cm

## ACCESSORIES



**Strap for pommel**  
862104-x  
For support over the knee

Two sizes



**Knee supports, fixed**  
861144-x  
Supports can be moulded to fit precisely without using tools

Size	Measurements
1	For size 1+2
2	For size 3+4



**Knee supports, hinged, fixed**  
861144-xx  
Provides support above and below the patella

Size	Measurements
1	For size 1+2
2	For size 3+4



**Knee supports, swingaway**  
861154-x Can be moved away from child without removing from the frame

Size	Measurements
1	For size 1+2
2	For size 3+4



**Knee supports, hinged, swing-away**  
861154-xx Provides support above and below the patella. Can be moved away from child without removing from the frame

Size	Measurements
1	For size 1+2
2	For size 3+4



**Formed foot plates, set**  
8614x  
For both prone and supine standing

Size	Measurements
1	L: 21½ cm (for sz. 1+2)
2	L: 24 cm (for sz. 3+4)



**Heel supports**  
86125x  
For both prone and supine standing

Size	Measurements
1	For size 1+2
2	For size 3+4



**Tray**  
862130-x Plastic tray with bowl and cover. Can be used for both Prone and Supine Standing

Size	Measurements
2	For size 1+2
3	For size 3+4



**Arm rests for plastic tray**  
861131  
Specially used for supine standing. One size fits all



**Toybar for plastic tray**  
861135  
Specially used for prone or supine standing. One size fits all



**Grip handle for plastic tray**  
861136  
Specially used for prone or supine standing. One size fits all