Nimbo Lightweight Posterior Posture Walker





- » Aluminium Frame in bright colours.
- » Height adjustable.
- Soft rubber wheels adhere to any surface and will not allow the Nimbo to slide backwards.
- » One directional rear wheels with aluminium ratchet and pin mechanism.
- » Anti-Reverse override bracket disengages the one directional rear wheels and allows forward & reverse mobility.
- » Front swivel wheels enhance maneuverability and ease of turning.
- » Locking mechanism easily switches front wheels from swivel to non-swivel.
- » Each wheel is independently adjusted.
- » Skid resistant, non-marking and easy rolling wheels.
- » Latex-free handgrips with flanged ends.
- » Folds easily for convenient transporting.
- » 2 year Limited Warranty.





Helps the user participate in outdoor activities. Ample 9" diameter height adjustable wheels for performance.

CODE	DESCRIPTION
KA1200N	Posture Walker - Tyke, Goldenrod Yellow
KA2200N	Posture Walker - Junior, Lime Green
KA3200N	Posture Walker - Youth, Cornflower Blue
KA4200N	Posture Walker - Young Adult, Midnight Blue
KA5200N	Posture Walker - Adult, Flame Red

CODE	HANDLEBAR HEIGHT RANGE	OVERALL WIDTH	WIDTH INSIDE HANDLE GRIP	DEPTH OPENED AT BASE	WEIGHT	WEIGHT CAPACITY
KA1200N	37-47cm (14½" - 18½")	53cm (21")	34cm (13½")	64cm (25")	4kg (9lb)	45kg (7st)
KA2200N	48-64cm (19" - 25")	56cm (22")	34cm (13½")	74cm (29")	5kg (11lb)	45kg (7st)
KA3200N	50-77cm (23" - 30½")	64cm (25")	38cm (15")	84cm (33")	5kg (12lb)	90kg (14st)
KA4200N	72-91cm (28½" - 36")	64cm (25½")	41cm (16")	91cm (36")	6kg (13lb)	90kg (14st)
KA5200N	92-104cm (36" - 41")	74cm (29")	48cm (19")	94cm (37")	7kg (16lb)	125kg (20st)



Nimbo Forearm Platform

Constructed of lightweight aluminium, the Nimbo forearm platforms are height and depth adjustable. The handgrips extend forward or back to the position needed, supinate, pronate and are angle adjustable. Each side is individually adjustable and allows custom positioning for each user.

The contoured armrest is moulded in one piece with a flexible edge for added comfort and has Velcro straps to prevent arms from slipping. They can be mounted to the mounting bracket either towards the inside or outside of the walker for lateral adjustment. Spacers are included for additional lateral range.



Pelvic Stabalizer

The pelvic stabilizer is adjustable in width/height and features back & lateral pads and a support belt. It is recommended for users who have an asymmetrical posture when standing or walking and therefore, cannot stay centered in the walker.

Seat

Seat mounts on standard walker. They fold up for standing and walking and fold down for convenience.



Adjustable sling seat supports the user that cannot stand completely on their own.



CODE	DESCRIPTION
KA1035FPS	Forearm Platform + Mounting Bracket - Small - 1pr/bx
KA1035FPL	Forearm Platform + Mounting Bracket - Large - 1pr/bx
KA1050	Pelvic Stabilizer - For use with all Nimbo Posture Walkers
KA1285	Seat - For use with KA1200N
KA2285	Seat - For use with KA2200N
KA3285	Seat - For use with KA3200N
KA4285	Seat - For use with KA4200N
KA5285	Seat - For use with KA5200N
CE1070S	Seat Harness - For use with KA 1200, KA 2200 and KA 3200
CE1070L	Seat Harness - For use with KA 4200
KA8100	Outdoor Wheels - For use with KA2200N - KA5200N
KA8100S	Outdoor Wheels - For use with KA1200N